



MODEL 8450
BABY SCALE
Operation Manual



8535-M151-O1 Rev A
10/02

CARDINAL SCALE MFG. CO.
PO BOX 151 • WEBB CITY, MO 64870
PH (417) 673-4631 • FAX (417) 673-5001
Web Site - <http://www.detectoscale.com>

Printed in USA

INTRODUCTION

Thank you for purchasing our Detecto Model 8450 Baby Scale. The 8450 comes complete with a one-piece platform and base with padded handles. With a 9V DC adapter (not included) the scale can be used in a fixed location or with a 9V battery (not included) as a portable scale.

This manual will guide you through installation, and operation of your scale. Please read it thoroughly before attempting to operate this scale and keep it handy for future reference.

FCC COMPLIANCE STATEMENT

WARNING! This equipment generates, uses and can radiate radio frequency and if not installed and used in accordance with the instruction manual, may cause interference to radio communications. It has been tested and found to comply with the limits for a Class A computing device pursuant to Subpart J of Part 15 of FCC rules, which are designed to provide reasonable protection against such interference when operated in a commercial environment. Operation of this equipment in a residential area may cause interference in which case the user will be responsible to take whatever measures necessary to correct the interference.

You may find the booklet "How to Identify and Resolve Radio TV Interference Problems" prepared by the Federal Communications Commission helpful. It is available from the U.S. Government Printing Office, Washington, D.C. 20402. Stock No. 001-000-00315-4.

All rights reserved. Reproduction or use, without expressed written permission, of editorial or pictorial content, in any manner, is prohibited. No patent liability is assumed with respect to the use of the information contained herein. While every precaution has been taken in the preparation of this manual, the Seller assumes no responsibility for errors or omissions. Neither is any liability assumed for damages resulting from use of the information contained herein. All instructions and diagrams have been checked for accuracy and ease of application; however, success and safety in working with tools depend to a great extent upon the individual accuracy, skill and caution. For this reason the Seller is not able to guarantee the result of any procedure contained herein. Nor can they assume responsibility for any damage to property or injury to persons occasioned from the procedures. Persons engaging the procedures do so entirely at their own risk.

Serial Number _____ Date of Purchase _____ Purchased Form _____
RETAIN THIS INFORMATION FOR FUTURE USE


TABLE OF CONTENTS

- Introduction 1
- Specifications 1
- Installation 2
 - Unpacking 2
 - Power Supply 2
 - Battery 2
 - Low Battery 3
- Display And Key Functions 3
- Features 3
- Operation 4
- Configuration 5
- Calibration 7
- Value Of Gravity For Calibration 8
- Display Messages 8


SPECIFICATIONS

Power Requirements	9V battery or 9V 100mA DC adapter <i>(not included)</i>
Capacity	44 lb x 0.02 lb, 704 oz x 0.5 oz (20 kg x 0.01 kg)
Weight Units	lb, oz, lb:oz, kg
Display	1.0" LCD display with 4 1/2 active digits
Operation Temperature	41 to 95 °F (5 to 35 °C)
Function Keys	ON/ZERO/OFF, Kg/Lb, MEMORY
Dimensions	Scale with handle: 24 x 14 x 2 in. (620 x 350 x 50 mm) Platform surface: 20 x 14 x 2 in. (500 x 350 x 50 mm)

PRECAUTIONS
 Before using this instrument, read this manual and pay special attention to all "WARNING" symbols:



IMPORTANT



ELECTRICAL WARNING

This is a “Table of Contents preview” for quality assurance

The full manual can be found at <http://the-checkout-tech.com/estore/catalog/>

We also offer free downloads, a free keyboard layout designer, cable diagrams, free help and support.

<http://the-checkout-tech.com> : *the biggest supplier of cash register and scale manuals on the net*